

# FEI CM30T TEM

## Instrument capabilities:

1. Instrument specifications:
  - a) Accelerating voltages: 100–300 kV
  - b) LaB<sub>6</sub> emitter
  - c) Resolution (at 300 kV): ~ 0.25 nm point; ~ 0.14 nm line
  - d) Minimum probe size: ~ 9 nm
2. Operating modes: CTEM, CBED, SAED, light element XEDS, PEELS, 1 Mpixel CCD camera (either 15 frames/second or slower single-image capture).
3. Specimen holders:
  - a) Double Tilt ( $\pm 60^\circ \alpha$ ,  $\pm 30^\circ \beta$ ):
    - with Be cup for XEDS
    - liquid N<sub>2</sub> cooled (96 K) with Be cup
    - heating (1270 K)
  - b) Tilt/rotate ( $\pm 60^\circ$ , rotation 360°)
  - c) Single Tilt ( $\pm 60^\circ$ )

## Typical experiments (examples):

- Quantitative XEDS
- Morphological and diffraction contrast studies of defects
- In situ heating & cooling studies
- Electron crystallography
- Weak beam studies of defects

